

|               |   |
|---------------|---|
| Product Name: | ABT-RFC-30M   |
|               |   |
| Description:  | <p>Standard Refrigerator &amp; Freezer Combination<br/> 30 cu. ft. capacity (25 ref./5 freezer)<br/> One exterior solid swing door and one interior solid swing door, right hinged<br/> Mechanical thermostat<br/> Refrigerator temperature range: 0°C to 10°C<br/> Freezer temperature range -15°C to -20°C<br/> 3 adjustable shelves, refrigerator<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> Stainless steel door handle<br/> Leveling legs<br/> Two probe access ports (3/8"); one, refrigerator/one, freezer<br/> Powder coat paint exterior<br/> High density urethane foam cabinet and door insulation<br/> CFC-free refrigerant (R404a), foam insulation and packaging<br/> Two year parts and labor warranty, plus an additional three year compressor parts warranty<br/> Exterior dimensions: 35 ¾ "W x 35 ½ "D x 72 ½ "H<br/> Manual defrost<br/> 115V, 60 Hz, 6 Amps, 1/3 HP<br/> Shipping weight: 383 lbs.</p> |