

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Research[®] plus Pipettes

Eppendorf Research [®] plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL	
epT.I.P.S. [®]	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention	■			■					■											
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■		■		■	■	■	■		
ep Dualfilter T.I.P.S. [®] LoRetention	■			■																
ep Dualfilter T.I.P.S. [®] SealMax	■			■				■	■	■	■									
Fixed volume, single channel pipettes	■																			
10 µL	■	■	■	■																■
20 µL	■	■	■	■																■
10 µL	■						■	■	■											
20 µL	■						■	■	■	■										
25 µL/50 µL/100 µL	■							■	■	■										
200 µL	■								■	■										
200/250/500/1,000 µL	■										■	■	■							
Variable volume, single-, 8- and 12-channel pipettes	■																			
0.1 – 2.5 µL	■	■	■	■																
0.5 – 10 µL	■	■	■	■																■
2 – 20 µL	■	◆	◆	■	■															■
2 – 20 µL	■						■	■	■	◆										
10 – 100 µL	■							◆	■	■	■									
20 – 200 µL	■								◆	◆	■									
30 – 300 µL	■									◆	◆	■								
100 – 1,000 µL	■										■	■	■							
120 – 1,200 µL	■											■	■							
0.25 – 2.5 mL	■													■						
0.5 – 5 mL	■														■	■				
1 – 10 mL	■																■	■		
Variable volume, 16- and 24-channel pipettes	■																			
1 – 20 µL	□				■															
5 – 100 µL	■					■														

■ normal volume ◆ limited volume



For more information go to www.eppendorf.com/eptips

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®] and Eppendorf Research[®] are registered trademarks of Eppendorf AG, Germany. GEloader™ and Microloader™ are trademarks of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Reference[®] 2 Pipettes

Eppendorf Reference [®] 2	0.1 – 10 µL S dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL
epT.I.P.S. [®]	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention	■			■					■								
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■				■	■	■		
ep Dualfilter T.I.P.S. [®] LoRetention	■			■					■								
ep Dualfilter T.I.P.S. [®] SealMax	■			■					■								
Fixed volume																	
1 µL/2 µL	■	■	■														■
5 µL/10 µL	■	■	■	■													■
20 µL	■		■	■													
10 µL	■				■	■	■										
20 µL	■				■	■	■	■									
25 µL/50 µL/100 µL	■					■	■	■									
200 µL	■						■	■									
200/250/500/1,000 µL	■								■	■	■						
2 mL/2.5 mL	■											■					
Adjustable volume																	
0.1 – 2.5 µL	■	■	■	■													
0.5 – 10 µL	■	■	■	■													■
2 – 20 µL	■	◆	◆	■													■
2 – 20 µL	■				■	■	■	◆									
10 – 100 µL	■				◆	■	■	■									
20 – 200 µL	■				◆	◆	■	■									
30 – 300 µL	■				◆	◆	◆	■									
100 – 1,000 µL	■								■	■	■						
0.25 – 2.5 mL	■											■					
0.5 – 5 mL	■												■	■			
1 – 10 mL	■														■	■	

■ normal volume ◆ limited volume



For more information go to www.eppendorf.com/eptips

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®] and Eppendorf Reference[®] are registered trademarks of Eppendorf AG, Germany. GEloader™ and Microloader™ are trademarks of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Xplorer[®]/Xplorer[®] plus Pipettes

Eppendorf Xplorer [®] / Eppendorf Xplorer [®] plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5 – 20 µL	
epT.I.P.S. [®]	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention	■			■																
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■		■		■	■	■	■		
ep Dualfilter T.I.P.S. [®] LoRetention	■			■																
ep Dualfilter T.I.P.S. [®] SealMax	■			■				■	■	■	■									
Variable volume, single-, 8- and 12-channel pipettes	■																			
0.5 – 10 µL	■	■	■	■																■
1 – 20 µL	■	◆	◆	■	■															
5 – 100 µL	■							◆	■	■										
10 – 200 µL	■							◆	◆	■										
15 – 300 µL	■							◆	◆	◆										
50 – 1,000 µL	■										■		■							
50 – 1,200 µL	■											■	■							
0.125 – 2.5 mL	■													■						
0.2 – 5 mL	■														■	■				
0.5 – 10 mL	■																■	■		
Variable volume, 16- and 24-channel pipettes	■																			
1 – 20 µL					■															
5 – 100 µL						■														

■ normal volume ◆ limited volume



For more information go to www.eppendorf.com/eptips

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®] and Eppendorf Xplorer[®] are registered trademarks of Eppendorf AG, Germany. GEloader™ and Microloader™ are trademarks of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.